

**A.R.C.  
LASER**

# FOX ENDO-PHOTOACOAGULATION

1.5 W

IDEAL IN  
CLINIC &  
OFFICE



Ophthalmology



### FOX 514 nm, 1.5 WATT

Amount	Order #	Description
1 x	LS03200- FAAB	FOX Diode Laser 514 nm red aiming beam foot switch 2




### ENDO COAGULATION use

#### CONSUMABLES C

<p>C PCS=5   LL13006s</p> <p>2</p> <p>Endo Probe, 23 G straight, sterile</p>	<p>C PCS=5   LL13025s</p> <p>3</p> <p>Endo Probe, 23 G curved, sterile</p>	<p>C PCS=5   LL13010s</p> <p>2</p> <p>Endo Probe, 25 G straight, sterile</p>
<p>C PCS=5   LL13015s</p> <p>3</p> <p>Endo Probe, 25 G curved, sterile</p>		

### EYE SAFETY

AS01003 safety goggle PCS=1



445/514/532/  
980/1064 nm



BG06076  
PCS=1  
514 nm,  
Zeiss



BG06077  
PCS=1  
514 nm,  
Leica

### EXCLUSIVE ACCESSORY



VP03104  
Trolley Ready Case



WE01447  
Trolley Cart